

## **AlloAid® DIP**

### **INDICATIONS**

AlloAid® DIP allograft is a Human Cellular and Tissue Based Product (HCT/P) per 21 CFR Part 1271. AlloAid® DIP sterile machined cortical bone allografts are intended for transplant in small bone fusion procedures. Each allograft is restricted to homologous use for transplant in fusion surgical procedures on a single occasion by a licensed physician or surgeon.

### **WARNINGS AND PRECAUTIONS**

An allograft may not elicit proper response from the recipient (e.g. fusion/union with adjacent tissue). It is possible for a host site to become infected. The allograft may not provide mechanical support and collapse, or may cause an inflammatory response. Current technologies may not preclude the transmission of infectious agents or disease, including hepatitis and HIV.

AlloAid® DIP may contain trace amounts of one or more processing agents including iodine, ethanol, hydrogen peroxide, gentamicin and/or vancomycin. It should not be used in patients with sensitivities to these processing agents.