

NeoSpan® SE

INDICATIONS

The In2Bones NeoSpan® Compression Staple Implant w/instruments is indicated for hand and foot bone fragments, osteotomy fixation and joint arthrodesis.

CONTRAINDICATIONS

General contraindications for the use of these implants for fusion, osteotomy and arthrodesis include:

- Significant bone demineralization
- Inadequate neurovascular status
- Inadequate skin or musculotendinous system
- Inadequate bone stock
- Physiologically or Psychologically unsuitable patient
- Possibility for conservative treatment
- Bone, musculature, tendons, or adjacent soft tissue compromised by disease, infection, or prior implantation, which cannot provide adequate support or fixation for the prosthesis
- Known metal allergy
- Diabetes
- Active infection
- Possibility of conservative treatment
- Growing patients with open epiphyses
- Patients with high level of activity