

Reference Toe System®

INDICATIONS

The RTS® Flexible 1ST MPJ Implant w/Grommets is intended for use in the treatment of:

- Hallux limitus or hallux rigidus
- Painful rheumatoid arthritis
- Hallux abducto valgus associated with arthritis
- Unstable or painful joint from previous surgery

CONTRAINDICATIONS

General contraindications for the use of these implants for joint reconstruction, osteotomy or fusion include:

- Significant bone demineralization.
- Inadequate neurovascular status.
- Inadequate skin or musculotendinous system.
- Inadequate bone stock.
- Physiologically or Psychologically unsuitable patient.
- Possibility for conservative treatment.
- Bone, musculature, tendons, or adjacent soft tissue compromised by disease, infection, or prior implantation, which cannot provide adequate support or fixation for the prosthesis.
- Known metal allergy. (1st MPJ w/GrommetsSystem Only)
- Diabetes.
- Active infection.
- Growing patients with open epiphyses.
- Patients with high level of activity