

Vitrium® BGS Wedge

INDICATIONS

Bone Graft Substitute:

The BGS implant is intended for use as a bone void filler for bony voids or gaps that are not intrinsic to the stability of the bony structure. BGS is indicated to be gently packed into bony voids or gaps of the skeletal system (i.e. extremities and pelvis). These defects may be surgically created osseous defects or osseous defects created from traumatic injury to the bone. The product provides a bone void filler that resorbs and is replaced with bone during the healing process.

CONTRAINDICATIONS

Bone Graft Substitute BGS is contraindicated:

- Possibility for conservative treatment
- When the device is intended as structural support in load-bearing bone
- Uncooperative patients who will not or cannot follow postoperative instructions, including individuals who abuse drugs and/or alcohol
- Inadequate skin, bone or neurovascular status
- Infection
- Irreparable tendon system
- Growing patients with open epiphyses
- Patients with high levels of activity
- It is not intended for use in the disc space