

i2b™ Rapid Set

INDICATIONS

i2b™ Rapid Set is an osteoconductive putty that is intended to be used to fill voids and gaps that are not intrinsic to the stability of the bone structure. These gaps or voids may be located in the extremities, spine, pelvis, or cranium. The putty may be shaped and pressed into the void by hand or inserted into a syringe and injected into the surgical site. The paste set in situ or ex situ provides a void filler that can augment hardware to support bone fragments during the surgical procedure. The set putty acts as a temporary support medium and is not intended to provide structural support during the healing process. The implant is radio-opaque. It is biocompatible and resorbs in the body as bone in-growth occurs.

CONTRAINDICATIONS

i2b™ Rapid Set is not designed or sold for any use except as indicated. Do not use i2b™ Rapid Set in the presence of any contraindication. The implant will not function if the implanted site is not well vascularized. The implant is contraindicated where it is intended as structural support in the skeletal system (e.g. mandibular segment replacement).