

## **Hercules® Suture Anchor System**

### **INDICATIONS**

The In2Bones Hercules® FT Suture Anchors are intended for use in the following applications:

1. Shoulder: Rotator Cuff, Bankart and SLAP lesion repair, Biceps tenodesis, Acromio-Clavicular separation and Deltoid repair, Capsular shift and Capsulolabral reconstruction.
2. Foot/Ankle: Lateral and Medial stabilization, Achilles tendon and Metatarsal ligament repair, Hallux Valgus and Midfoot reconstruction.
3. Knee: Medial collateral and Lateral collateral ligament repair, Patellar tendon and Posterior oblique ligament repair.
4. Hand/Wrist: Scapholunate ligament, Radial collateral ligament and Ulnar collateral ligament reconstruction.
5. Elbow: Biceps tendon reattachment, Tennis elbow repair, Ulnar and Radial collateral ligament reconstruction.

### **CONTRAINDICATIONS**

General contraindications for the use of these implants for joint reconstruction include:

- Significant bone demineralization.
- Inadequate neurovascular status.
- Poor blood supply, which may impair the healing process
- Inadequate skin or musculotendinous system.
- Inadequate bone stock.
- Physiologically or Psychologically unsuitable patient.
- Surgical procedures other than those indicated.
- Possibility for conservative treatment.
- Bone, musculature, tendons, or adjacent soft tissue compromised by disease, infection, or prior implantation, which cannot provide adequate support or fixation for the prosthesis.
- Known allergies.
- Diabetes.
- Active infection to materials.
- Patients with high level of activity.
- Pathological conditions of the soft tissue which would compromise secure fixation.