

INDICATIONS

AlloAid® SubChon Implant is composed of 100% human bone tissue and produced from segmented cancellous bone.

The graft is indicated for use in bony voids or gaps that are not intrinsic to the stability of the bony structure. It should be gently packed into bony voids or gaps of the skeletal system. These defects may be surgically created osseous defects or osseous defects created from traumatic injury to the bone. The graft provides a bone void filler that resorbs and is replaced by the growth of new bone during the healing process. It can be mixed with autogenous bone marrow prior to use at the physician's discretion.

CONTRAINDICATIONS

Do not use the graft in the presence of any contraindication. The graft is contraindicated where the graft is intended as structural support in the skeletal system. Other conditions representing relative contraindication include:

- Severe vascular or neurological disease
- Uncontrolled diabetes
- Pregnancy
- Uncooperative patients who cannot or will not follow post-operative instruction, including individuals who abuse drugs and/or alcohol
- Hypercalcemia, abnormal calcium metabolism
- Existing acute or chronic infections, especially at the site of the operation
- Inflammatory bone disease such as osteomyelitis
- Malignant tumors
- Severely impaired renal function